

Immunofluorescence Staining of Cells for Flow Cytometry

Overview

The following is a general cell staining protocol for analysis by flow cytometry. It is recommended that a separate isotype control be tested in parallel as a negative control.

Material Required

- Flow Buffer
- Trypan blue
- 10% formalin solution
- PBS
- 15 mL conical tubes
- 12 × 75 mm round-bottom tubes

Material Preparation

Flow Buffer

- 99 mL PBS (1×)
- 1 mL FBS
- 0.1 mL Sodium Azide (100%)
 - If detecting **extracellular antigens**, the Flow Buffer is ready for use.
 - If detecting **intracellular antigens**, add Saponin to a final concentration of 0.1%.
- Store Flow Buffer at 4°C.

Stain Cells for Flow Cytometry

1. Isolate and dissociate the cells to a single cell suspension. Collect the cells in a 15 mL conical tube.
2. Remove a sample of the cell suspension and count live cells using trypan blue and a hemocytometer.
3. Centrifuge the cell suspension at $300 \times g$ for 5 minutes at 4°C .
4. Aspirate the medium and flick the tube to disrupt the cell pellet.
5. If detecting **extracellular antigens**, resuspend the cell pellet in 3 mL of PBS.
6. If detecting **intracellular antigens**, resuspend the cell pellet in 2 mL 10% formalin solution and incubate for 15 minutes at room temperature.

Note: Formalin should be added slowly while agitating the tube to avoid cell clumping.

7. Centrifuge the cell suspension at $300 \times g$ for 5 minutes at 4°C .
8. Aspirate the medium and flick the tube to disrupt the cell pellet.
9. Add enough Flow Buffer to bring the cell suspension to a concentration of 2×10^6 to 1×10^7 cells per mL. Keep the cells on ice.

Note: The Flow Buffer formulation depends on if the antigen is extra or intracellular. Refer to the preparation section to make sure the appropriate Flow Buffer is used.

10. For each sample, add 100 μL of the cell suspension to a 12×75 mm round-bottom tube.
11. Add the appropriate amount of primary antibody and/isotype control to each sample.
12. Incubate on ice for 30 minutes to 1 hour, protecting the samples from light if using a conjugated antibody.
13. Add 4 mL of Flow Buffer.
14. Centrifuge at $300 \times g$ for 5 minutes at 4°C .
15. Aspirate the supernatant and flick the tube to disrupt the cell pellet.

Note: if using a conjugated antibody, skip steps 16 through 20 and go directly to step 21. If using a purified antibody, continue with step 16.

16. Add the appropriate amount of secondary conjugated antibody to each sample.
17. Incubate on ice for 30 minutes to 1 hour, protecting the samples from light.
18. Add 4 mL of Flow Buffer
19. Centrifuge at $300 \times g$ for 5 minutes at 4°C .
20. Aspirate the supernatant and flick the tube to disrupt the cell pellet.
21. Add appropriate volume of Flow Buffer to each tube.

*Note: For this step, use Flow Buffer **without** 0.1% Saponin for both extra and intracellular antigen detection.*

22. Analyze the cells by flow cytometry within 4 hours.

REPROCELL USA Inc

9000 Virginia Manor Road
Suite 207
Beltsville, MD 20705
USA

Tel: +1 301 470 3362

Email: info-us@reprocell.com

REPROCELL Europe Ltd

Thomson Pavilion, Todd Campus
West of Scotland Science Park
Acre Road Glasgow
G20 0XA
UK

Tel: +44 (0)141 465 3460

Email: info-emea@reprocell.com

REPROCELL India Ltd

3-1-135/1A, CNR Complex
Mallapur
Hyderabad 500 076
Telangana
India

Tel: +44 (0)141 465 3460

Email: Bhargavi.Gurram@reprocell.com

REPROCELL Inc (Japan)

MetLife Shin-yokohama 381, Bldg. 9F
3-8-11, Shin-yokohama Kohoku-ku,
Yokohama
Kanagawa 222-0033
Japan

Tel: +81 45 475 3887

Email: info-asia@reprocell.com