

## Recombinant Human Laminin Fragments for Cell Culture







### Overview

Culturing human stem cells and iPS cells under feed-free conditions requires the use of extracellular matrix proteins (ECM) as an anchor to promote adherence to laboratory plastic-ware. Laminin fulfills this purpose, and through its binding of membrane bound integrin, multiple intracellular signal cascades are activated. The recombinant human laminin-511 E8 fragment has become the most popular ECM for human stem cell research due to its broad applicability and strong capacity to induce cell proliferation. iMatrix-211 is a widely used matrix for culture of cardiomyocytes.



### Matrixome® iMatrix™ Substrates available from REPROCELL

- Recombinant human protein produced in CHO-S cells with serum-free medium
- iMatrix-511 SILK produced in Silkworms
- Ready-to-use liquid (frozen) format
- Each lot validated for high performance in human cell culture
- Stamped with an expiration date
- Tested for endotoxin, mycoplasma and bacterial contamination
- Guaranteed more than 95% pure Laminin
- Integrin binding activity quality checked

Product	Description	Typical Usage	Cat. No.	Size
 iMatrix-221	Recombinant Human Laminin-211 E8 Fragments	Cardiomyocytes	NP892-061 NP892-062	2 × 175 µg 6 × 175 µg
 iMatrix-411	Recombinant Human Laminin-411 E8 Fragments	Endothelial cells	NP892-041 NP892-042	2 × 175 µg 6 × 175 µg
 iMatrix-511	Recombinant Human Laminin-511 E8 Fragments	Pluripotent Stem Cells	NP892-011 NP892-012	2 × 175 µg 6 × 175 µg
 iMatrix-511 Silk	Recombinant Human Laminin-511 E8 Fragments, expressed in Silkworm	Pluripotent Stem Cells	NP892-021	2 × 175 µg

## What is Laminin E8?

Laminin is a ubiquitous cell surface protein comprised of multiple sub-units. E8 fragments are proteolytic fragments that retain the high binding capacity of full-length laminin. At least 15 different sub-types of laminin have been discovered. For example, laminin-511 is comprised of the  $\alpha 5$ -chain,  $\beta 1$ -chain, and  $\gamma 1$ -chain. This sub-type is known to bind strongly to  $\alpha 6\beta 1$  integrin.



## iPS Cell Colony Morphology

Human iPS cell lines 201B7 (retrovirus reprogrammed) and RC010 (mRNA reprogrammed) are shown grown on iMatrix-511 coated plates in StemFit™ AK02 medium. Both exhibit flat, rounded colonies with distinct edges; characteristic of healthy, pluripotent iPS cell colonies.

## NEW iPS Cell Reprogramming System

The NEW StemRNA™-3<sup>rd</sup> Gen Reprogramming Kit from Stemgent® (00-0076) used in combination with iMatrix-511 and NutriStem™ hPSC XF medium, brings you higher efficiencies, more options, and simple, improved protocols. One kit enables you to reprogram primary human fibroblasts, blood-derived endothelial progenitor cells (EPCs), or urine-derived epithelial cells (UDCs). The entire workflow is xeno-free with cGMP compatible reagents. Delivered as a mRNA cocktail, the reprogramming factors enable the high efficiency conversion to iPS cells in as little as 2 weeks without the risk for vector-induced genome mutations. RNA-based approaches are regarded as one of the most promising clinically compatible methods for iPS cell reprogramming.



Cat. No.	Description	Size
00-0076	StemRNA™ 3rd Gen Reprogramming Kit	1 kit
01-0005	NutriStem hPSC XF Culture Medium	500 mL
01-0020-50	NutriFreez™ D10 Cryopreservation Medium	50 mL
ASB01	StemFit™ Basic02 Medium	500 mL
04-0012	Stemolecule™ Y27632	2 mg

## References

1. Ido H. et.al. *J. Biol. Chem.* 282(15): 11144-54, 2007
2. Taniguchi Y. et.al. *J. Biol. Chem.* 284(12) 7820-31, 2009
3. Miyazaki T. et.al. *Nature Commun.* 3: 1236, 2012
4. Nakagawa M. et.al. *Sci. Rep.* 4: 3594, 2014
5. Doi D. et.al. *Stem Cell Reports* 2(3): 337-50, 2014
6. Takashima Y. et.al. *Cell* 158(6): 1254-69, 2014
7. Fukuta M. et.al. *PLoS One* 9(12): e112291, 2014

Matrixome company name and logo and iMatrix brand name are the property of Matrixome Corp., Japan.  
Company names, brands, logos and trademarks for REPROCELL, BioServe, Stemgent, Alvetex and Biopta are all the property of REPROCELL Inc.  
© 2019 REPROCELL, Inc. All rights reserved.

### REPROCELL Europe Ltd

Thomson Pavillion, Todd Campus  
West of Scotland Science Park  
Acre Road  
Glasgow, G20 0XA  
UK

T: +44 (0)141 465 3460  
E: info-emea@reprocell.com

### REPROCELL USA Inc

9000 Virginia Manor Road  
Suite 207  
Beltsville, MD 20705  
USA

T: +1 301 470 3362  
E: info-us@reprocell.com

### REPROCELL India Ltd

3-1-135/1A, CNR Complex  
Mallapur  
Hyderabad 500 076  
Telangana  
India

T: +91-40-27178178  
E: Bhargavi.Gurram@reprocell.com

### REPROCELL Inc (Japan)

MetLife Shin-yokohama 381, Bldg. 9F  
3-8-11, Shin-yokohama  
Kohoku-ku, Yokohama  
Kanagawa 222-0033  
Japan

T: +81 45 475 3887  
E: info-asia@reprocell.com