

**RESEARCH USE ONLY**

Ver. 001

Product No. 892 441: 1 vial

Product No. 892 442: 2 vials

**Purified recombinant thermolysin  
Brightase-TH (4.7 mg/vial)****Product Description:**

Brightase-TH, the recombinant thermolysin from *Bacillus thermoproteolyticus rokko*<sup>1)</sup>, is produced with stable activity using the *Brevibacillus* expression system<sup>2)</sup>. This system leads to the expression of recombinant proteins with low endotoxin contamination since *Brevibacillus* is a Gram-positive bacterium. Our purified thermolysin has typical neutral proteinase activity. Due to no collagenolytic activity in Brightase-TH, it is necessary to add Brightase-C (recombinant collagenase) for dissociation of tissues to isolate primary cells.

**Content:**

Recombinant thermolysin from *Bacillus thermoproteolyticus rokko* (lyophilized): 4.7 mg/ vial

- ※ Lyophilized at 0.94 mg/ml of 5 mL solution  
(HEPES buffer: 50 mM HEPES, 2 M NaCl, 5 mM CaCl<sub>2</sub>, pH 7.5)
- ※ Organism of origin: *Bacillus thermoproteolyticus rokko*<sup>1)</sup>
- ※ Recombinant production in: *Brevibacillus choshinensis*<sup>2)</sup>  
(Endotoxin-free ( $\leq 10$  EU/mg), animal-derived components free)

**Specification**

EC number	: 3.4.24.27
Activity*	: 0.8~1.2 MOCac U/ mg
Optimal pH	: 7.0~7.5
Recommended temperature for use:	30~40 °C
Thermal stability	: 80% of activity remains after enzyme is incubated in HEPES buffer at 70°C for 2 hours.
Inhibitors	: EDTA, phosphoramidon
M.W.	: 34,600
E <sup>1%</sup> <sub>1cm</sub> (280nm)	: 17.65 (pH7.0)
Content	: 4.7 mg/ vial
Endotoxin	: $\leq 10$ EU/ mg

\*One unit of activity was defined as the amount of enzyme releasing 1  $\mu$ M of MOCac-PLG from synthetic substrate (MOCac-PLGL(Dpa)AR) at pH 7.5 at 30°C per sec. <sup>3,4)</sup>

**Storage condition**

Brightase-TH should be stored at under -20 °C and is stable for 2 years at this condition. After dissolving a lyophilized product, the solution should be stored in aliquots at under -20 °C. Please avoid multiple freeze-thaw cycles.

**Stock solution**

Stock solution of Brightase-TH can be prepared in 5 mL sterile distilled water by agitating at 4°C for 30 min. It contains 0.94 mg/ml enzyme, 50 mM HEPES, 2 M NaCl and 5 mM CaCl<sub>2</sub> (pH 7.5) after reconstitution.

**Working solution**

The working solution is prepared by mixing the cold enzyme stock solution with cold buffer (containing  $\geq 1$  mM CaCl<sub>2</sub>). Calcium is necessary for activity of thermolysin. If sterile filtration is desired, a sterile filter with low protein binding properties should be used.

### Measurement of enzyme activity

Brightase-TH has a MOCAC-PLGL(Dpa)AR degrading activity, 0.8~1.2 MOCAC U/ mg at pH7.5 at 30°C per sec<sup>3, 4</sup>. One MOCAC unit is comparable to 100 caseinolytic activity (CU). One CU will hydrolyze casein to produce color equivalent to 1.0 μmol of tyrosine at pH7.5 at 37°C per min (Color by Folin-Ciocalteu reagent)<sup>5</sup>).

### References

- 1) O'Donohue, MJ. et al., *Biochemical Journal*, **300**, 599–603 (1994)
- 2) Mizukami, M. et al., *Current Pharmaceutical Biotechnology*, **11**, 251-258 (2010)
- 3) Muta, Y. and Inouye, K., *Journal of Biochemistry*, **132**, 945-951 (2002)
- 4) Oneda, H. et al., *Biosci. Biotechnol. Biochem.*, **68**, 1811-1813 (2004)
- 5) Anson, ML., *J. Gen Physiol*, **22**, 79-89 (1938)

### Contact and Support

E-mail: [protein-info@nippi-inc.co.jp](mailto:protein-info@nippi-inc.co.jp)

### Manufacturer

Nippi, Incorporated

Biological and Chemical Products Division

Head Quarter: 1-1-1 Senjumidori-cho, Adachi-ku, Tokyo JAPAN 120-8601

URL: <http://www.nippi-inc.co.jp/>