

KOKEN

AteloCell[®]

DATA SHEET

Atelocollagen coated β -TCP scaffold

Cat No. ACB-05S

Lot No.

DISC TYPE:	Honeycomb	Flat
DIAMETER :	3.0±0.3 mm	3.0±0.3 mm
HEIGHT :	1.0±0.1 mm	1.0±0.1 mm
HOLE: NUMBER:	37 holes	—
CROSS-SECTIONAL AREA:	0.07+0.03/-0.02 mm ²	—

	<u>Specification</u>	<u>Results</u>	<u>Results</u>
ATELOCOLLAGEN COAT: (Type I atelocollagen derived from Bovine dermis)	Coated	Pass	Pass
BIOLOGICAL INDICATOR TEST:	Negative	Pass	Pass

FOR RESEARCH USE ONLY, NOT FOR HUMAN BODY.

KOKEN Co., Ltd.
1-4-14 Koraku, Bunkyo-ku, Tokyo 112-0004, Japan
TEL: +81-3-3816-3525 FAX: +81-3-3816-3570

www.atelocollagen.com

ReproCELL, Inc.

KDX Shin-yokohama 381 Bldg 9F
3-8-11, Shin-yokohama, Kohoku ward,
Yokohama, Kanagawa 222-0033, Japan
Tel: +81 (0) 45 475 3887
Fax: +81 (0) 45 474 1006
Email: info_jp@reprocell.com