

**KOKEN**

AteloCell®

**DATA SHEET**

**Atelocollagen sponge, MIGHTY**

**Cat No. CSM-25**

**Lot No.**

ORIGIN:	Bovine dermis	
REFERENCES:	Refer to the AteloCell® website	
	<u>Specification</u>	<u>Results</u>
COMPRESSIVE STRENGTH:	10~30 kPa	19.76 kPa
SIZE		
DIAMETER:	4.3~4.7 mm	4.47 mm
HIGHT:	2.6~3.4 mm	2.98 mm
<p>Notes: These measurements are under dry condition. This product can be expanded to approximately 5 mm (diameter) and 3 mm (height) under wet condition.</p>		
STERILITY TEST: (Medium: TGC-I and SCD)	Negative	Pass
CELL CULTURE TEST: (Cell: Human Fibroblast)	Normal	Pass

FOR RESEARCH USE ONLY, NOT FOR HUMAN BODY.

KOKEN Co., Ltd.  
1-4-14 Koraku, Bunkyo-ku, Tokyo 112-0004, Japan  
TEL: +81-3-3816-3525 FAX: +81-3-3816-3570

**www.atelocollagen.com**

**ReproCELL, Inc.**

KDX Shin-yokohama 381 Bldg 9F  
3-8-11, Shin-yokohama, Kohoku ward,  
Yokohama, Kanagawa 222-0033, Japan  
Tel: +81 (0) 45 475 3887  
Fax: +81 (0) 45 474 1006  
Email: info\_jp@reprocell.com