



AteloCell[®]

DATA SHEET

Atelocollagen sponge, MIGHTY

Cat No. CSM-50

Lot No.

ORIGIN:	Bovine dermis	
REFERENCES:	Refer to the AteloCell [®] website	
	<u>Specification</u>	<u>Results</u>
COMPRESSIVE STRENGTH:	10~30 kPa	18.7 kPa
SIZE		
DIAMETER:	4.3~4.7 mm	4.43 mm
HEIGHT:	2.6~3.4 mm	2.78 mm
	(Notes : These measurements are under dry condition. This product can be expanded to approximately 5 mm (diameter) and 3 mm (height) under wet condition.)	
STERILITY TEST: (Medium: TGC-I and SCD)	Negative	Pass
CELL CULTURE TEST: (Cell: Mouse Chondrogenic Cell Line)	Normal	Pass

FOR RESEARCH USE ONLY, NOT FOR HUMAN BODY.

KOKEN Co., Ltd.

1-4-14 Koraku, Bunkyo-ku, Tokyo 112-0004, Japan
TEL: +81-3-3816-3525 FAX: +81-3-3816-3570

www.atelocollagen.com

ReproCELL, Inc.

KDX Shin-yokohama 381 Bldg 9F
3-8-11, Shin-yokohama, Kohoku ward,
Yokohama, Kanagawa 222-0033, Japan
Tel: +81 (0) 45 475 3887
Fax: +81 (0) 45 474 1006
Email: info_jp@reprocell.com